

Turn Your Existing GC–MS into an Online Organic Aerosol Analyzer

A modular front end for automated molecular analysis of organic aerosols.



Soft Connection

10 cm Needle +
Double Septa

No GC Modification

Easy Installation
& Switching

Automated Analysis

1–2 h Time Resolution
Continuous
Measurements

Broad Organic Coverage

Hydrocarbons +
Oxygenated Organics

EXISTING GC–MS



ACTS FRONT END



MOLECULAR-LEVEL GC–MS OUTPUT

How ACTS Works

Automated aerosol collection, thermal desorption and direct GC–MS introduction

MODE 1

Sampling + Analysis

Collection and Thermal Desorption Cell (CTD) · parallel sampling and GC–MS analysis

AEROSOL INLET → CYCLONE + HUMIDIFIER → CTD —
COLLECTION

NEEDLE BETWEEN SEPTA



GC–MS — PREVIOUS SAMPLE

- Aerosols are collected in the CTD
- Needle rests between the double septa
- Sampling and analysis proceed together

SAMPLING AND ANALYSIS RUN IN PARALLEL

MODE 2

Thermal Desorption

Collected organics are introduced directly into the GC–MS

He / MFC → CTD — THERMAL DESORPTION → 10 cm
NEEDLE

DOUBLE SEPTA + GC INLET



GC–MS

- Needle enters the GC inlet
- Collected organics are thermally desorbed
- Helium transfers analytes directly

10 cm DIRECT TRANSFER

Soft connection: 10 cm needle + double septa



- 01** Needle rests between the two septa during sampling.
- 02** Needle enters the GC inlet during thermal desorption.
- 03** No drilling, rigid metal unions or structural GC modification.

RAPID CONNECTION, REMOVAL AND SWITCHING

Designed for Existing GC–MS Platforms

	ACTS	Representative TAG-type configurations
CTD–GC connection	10 cm needle + double septa	High-temperature valve / rigid plumbing commonly used
GC modification	Not required	Dedicated integration may require modification
GC flexibility	Rapid connection / removal	More integrated dedicated configuration

TAG architectures have evolved; comparison refers to representative published configurations and illustrates system integration only.

Performance & Chemical Coverage

From n-Alkanes to Oxygenated Organics

ACTS Quantifies Oxygenated Organics Beyond Alkanes

~0.07 μm

Aerosol Collection D_{50}

>90%

Collection Efficiency
 $\geq 0.16 \mu\text{m}$

33–154%

C_{13} – C_{40} n-Alkane
Relative Transfer
Efficiency


~7%

7-Day Signal Drift

0.5–0.8 ng

C_{36} Hexatriacontane
Empirical Detection Limit

High Transfer Efficiency for Mid- to Low-Volatility Organics




Open squares CTD 350 °C · open circles CTD 320 °C · filled diamonds alkanoids acids, CTD 330 °C.
y-axis: relative transfer efficiency, CTD ÷ GC direct injection (%). x-axis: carbon number.


C_{36} Hexatriacontane = 154% at 350 °C

Measured CTD/direct-injection response; values >100% reflect the empirical relative-transfer method.


Beyond Alkanes: Oxygenated Organics




Decanoic acid · $R^2 \approx 0.99$



Undecanoic acid · $R^2 \approx 0.98$



Octadecanoic acid · $R^2 \approx 0.97$



Levoglucosan · $R^2 \approx 0.91$

Empirical Detection Limits

Compound	ACTS
Decanoic acid	$\leq 0.5 \text{ ng}$
Undecanoic acid	$\leq 0.5 \text{ ng}$
Octadecanoic acid	1–5 ng
Levoglucosan	5–10 ng

Sensitivity Benchmark

Compound	ACTS	Representative Published TAG
C_{16} Hexadecane	$\leq 0.05 \text{ ng}$	0.074 ng ¹
C_{36} Hexatriacontane	0.5–0.8 ng	7.2 ng ² · 0.5 ng ³
Decanoic acid	$\leq 0.5 \text{ ng}$	0.9 ng ¹
Levoglucosan	5–10 ng	6.3 ng ¹

ACTS measures both hydrocarbons and GC–MS-compatible oxygenated organics.

Comparable sensitivity to established TAG systems, with strong performance for selected high-carbon-number compounds.

Highly oxygenated species may require derivatization or other analytical approaches. TAG values represent selected published configurations and are not intended to represent all TAG systems. ¹ Williams et al., 2006; ² Kreisberg et al., 2009 (TAG–GC/MS); ³ Worton et al., 2012 (TAG–GC×GC/MS).

Applications in Real Aerosols


Molecular speciation across ambient, indoor and field aerosol samples

Styrene & Styrene Oligomers in Ambient Air

Shanghai Dianshan Lake Science Observation Station · 18 May 2023, 11:00

Styrene Monomer · Dimer · Trimer

AUTHENTIC RESOLVED CHROMATOGRAM




Numerous resolved molecular peaks demonstrate broad compound-level separation in a real ambient air sample.

Indoor Oxygenated Organic Aerosols

m/z 60 Extracted Ion Chromatograms · Morning · Noon · Afternoon


Shanghai Bay Valley Tech Park · September 2024



m/z 60 extraction highlights time-resolved changes in levoglucosan, alkanolic acids and other oxygenated organics.

Heavy Petroleum / UCM in Ambient Air

AMBIENT AIR · SAMPLE NO. 189 · 20 JUN 2026 · 22:00



Field evidence for heavy, low-volatility petroleum-related organics, with resolved peaks above the UCM envelope.

Molecular Speciation Across Sources and Functionalities

Source-specific tracers · Oxygenated organics · Heavy petroleum / UCM signatures

Peer-Reviewed Instrument Study

Yaping Zhang, Feng Shen, Yao Zhang, Lin Wang, Juntao Huo, Wenyu Chu, Xin Yang & Xiaofei Wang.
A low-cost in-situ instrument for organic aerosol speciation using thermal desorption. Aerosol Science and Technology (2026).
DOI 10.1080/02786826.2026.2614032

Yaping Zhang · EIPROTECH

contact@eiprotech.com
Phone / WeChat: +86 153 1779 0962
eiprotech.com

Explore ACTS



Specifications · Publications · Product Information · Contact